

Sperm CryoProtect™

Symbols

EN

Intended Use
Optimized for freezing human sperm.

Components	
Endotoxin free	EDTA
Potassium chloride	HEPES
Magnesium sulphate	Glucose
Potassium dihydrogen phosphate	Glycerol
Sodium bicarbonate	Purified water
Sodium pyruvate	
Calcium lactate	

Performance Characteristics	
pH	7,5-8,5
Endotoxin levels	<1,0 EU/ ml
Recovery rate of the original motile spermatozoa after freezing and thawing	>50%
Sterile filtered by human sperm survival only	SAL 10 ³
Contents are tested by human sperm survival only	
Bottles and caps are M.E.A. tested	

Storage and Stability
Store at 2 to 30°C and avoid temperatures above or below these values. Under these conditions Sperm CryoProtect has a shelf-life of 12 months. The expiry date is shown on both bottles and cartons.

Open and close bottles under aseptic conditions. After opening store at 2 to 8°C when not in use. Shelf-life on the product label applies when the product is stored and handled according to manufacturer's recommendations.

No antibiotics, unstable additives or preservatives have been added by the manufacturer to Sperm CryoProtect.

Precautions and Warnings

- Use aseptic procedures at all times
- Sperm CryoProtect contains glycerol which is combustible. A safety data sheet is available from the distributor or manufacturer (see www.nidacoin.com)
- Do not use any solution which shows evidence of bacterial contamination or if stopper accidentally comes in contact with unwashed surfaces
- Please check for regulatory compliance governing the use of ART products in your country
- Do not use contents if lamp-ember evident is broken
- Do not re-use
- Federal Law (USA) restricts this device to sale by or on the order of a physician

Please check for regulatory compliance governing the use of ART products in your country

Recommendations
Although it is possible to freeze unprocessed semen, we recommend that the ejaculate be prepared on a PureSperm density gradient and washed with PureSperm Wash before adding Sperm CryoProtect. This method removes seminal plasma as well as ROS and their sources, thereby ensuring optimal recovery of motile sperm after thawing.

Reagents and Equipment
• SpermCryoProtect and PureSperm Wash
• CryoFluor™
• Sterile pipettes
• Disposable sterile conical centrifuge tubes (e.g. Falcon 2075)
• Disposable sterile cryopreservation vials or plastic straws

Freezing sperm with Sperm CryoProtect

- Gradient-prepared sperm**
 - Add 1 part of Sperm CryoProtect to 3 parts of sample (e.g 100µL PureSperm CryoProtect to 300µL sperm sample),ensuring thorough mixing after adding each drop in order to avoid osmotic shock
 - Fill straws with equal suspension and seal the straws
 - Equilibrate straws in refrigerator for 10 - 60 minutes
 - Place the straws horizontally in nitrogen vapour on the CryoFluor™Straws.
- Close lid/ Leave for 10-30 minutes
- Transfer the straws quickly into the liquid nitrogen and store in liquid nitrogen. Do not touch the straws with your hand

- Thawing of gradient-prepared sperm**
 - Remove straw from liquid nitrogen tank
 - Place straw in water at 37°C for 30 secs
 - Dry surface of straw with sterile tissue paper
 - Cut one end of straw
 - Resuspend contents in 5ml PureSperm Wash by cutting the upper end of the straw. Any sperm suspension remaining in the straw can be expelled using a pipette
 - Centrifuge at 500 x g for 10 minutes. Do not use the brake
 - Aspirate PureSperm Wash supernatant leaving as much liquid as required for desired concentration. No pellet is seen, leave the bottom 100µL fluid

- Unprocessed ejaculate:**
 - Add 1 part of Sperm CryoProtect to 3 parts of sample (e.g 100µL PureSperm CryoProtect to 300µL sperm sample),ensuring thorough mixing after adding each drop in order to avoid osmotic shock
 - Transfer 0,8-1,8 ml of the mixture to 2 ml cryovials
 - Place the vials in the fridge (4-5°C) for 30 min
 - Freeze the vials horizontally in the freezer or in nitrogen vapour, above the liquid nitrogen surface on the CryoFluor™Vials. Close lid. Leave for 30 min
 - Store in liquid nitrogen

- Thawing of unprocessed ejaculate**
 - Remove the vials from liquid nitrogen tank
 - Place vials in water at 37°C until no ice crystals can be seen, approximately 2-3 min
 - Divide the thawed material with 0,5 ml PureSperm Wash
 - Prepare a 40% and 20% gradient with 1 ml, in each layer. Layer the thawed ejaculate onto the gradient
 - Centrifuge at 300 x g for 20 min
 - Aspirate everything except the pellet and add 4-6 ml of the PureSperm 80% layer
 - Use a new pipette to aspirate the pellet. Transfer to a new tube containing 4 ml PureSperm Wash
 - Resuspend contents in 5ml PureSperm Wash by cutting the upper end of the straw. Any sperm suspension remaining in the straw can be expelled using a pipette
 - Centrifuge at 500 x g for 10 min
 - Aspirate PureSperm Wash supernatant and the sample is now ready for use

To achieve the correct g force:
www.nidacoin.com/gp

Symbol	
1. Temperature limit	5. Consult instructions for use
2. Use by - see label	6. CE Mark
3. Sterilized using aseptic processing techniques	7. Manufacturer
4. Batch code	8. Intended Use

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